


# ETN专家级胫骨髓内钉

ETN Multi-direction Tibial Lockable Intramedullary Nail



- 多种治疗解决方案
- 一个产品满足胫骨近端、远端及远端股骨干骨折质量
- 更强的稳定性

- Several solutions for treatment
- xpert tibial nail is indicated for fractures of distal proximal and shaft.
- Stronger stability



# ETN专家级胫骨髓内钉

## ETN Multi-direction Tibial Lockable Intramedullary Nail

### 产品特点 Product features

- 先进的近端交锁设计, 为近端骨折块提供理想的稳定性  
 Advanced interlocking design at proximal side provides an ideal stability for all fragmented bones.
- 先进的远端交锁设计, 可应用于更远端的骨折  
 Advanced interlocking design at distal side can be applied to more distal fractures.
- 可通过尾帽固定一枚斜行锁钉来达到绝对的成角稳定  
 The end cap fix one oblique locking screw can achieve octet stability.
- 全部锁钉采用双线设计, 即使在松质骨中也能获得很好的把持力  
 The locking screws designed in double lines, even in cancellous bone also can get a good gripping force.

### 产品信息 Product Information

产品编号/Product No.	产品规格/Specification		材质/Material
	直径/Diameter	长度/Length(mm)	
4217.0825	8	255	钛合金 Titanium Alloy
4217.0827	8	270	
4217.0828	8	285	
4217.0830	8	300	
4217.0831	8	315	
4217.0833	8	330	
4217.0834	8	345	
4217.0836	8	360	
4217.0837	8	375	
4217.0925	9	255	
4217.0927	9	270	
4217.0928	9	285	
4217.0930	9	300	
4217.0931	9	315	
4217.0933	9	330	
4217.0934	9	345	
4217.0936	9	360	
4217.0937	9	375	
4217.1025	10	255	
4217.1027	10	270	
4217.1028	10	285	
4217.1030	10	300	
4217.1031	10	315	
4217.1033	10	330	
4217.1034	10	345	
4217.1036	10	360	
4217.1037	10	375	

### 适应症 Indications



### 苏州苏南捷迈得医疗器械有限公司

Suzhou Sunan Zimmered Medical Instrument Co.,Ltd.

地址: 江苏省张家港市合兴镇昆仑山路

Add: Kunlunshan road, Hexing, town, Zhangjiagang city, Jiangsu, province, China

国内 销售部电话: +86-512 5850 0520

China: +86-186 5112 8828

国际 Tel: +86-512 5691 3892

International: Tel/Whatsapp: +86-186 5112 8828



Wechat

传真Fax: +86-512-58508580

邮编(P.C): 215626

邮箱E-mail: jenny.wang@zimedede.com

网址Web: www.zimedede.com

企业微信号(Company WeChat ID): JMDYX1990

